# 3

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HOW TO EXERCISE  
WITH BOLTEX INERTIAL

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# How to Use

**Boltex**  
INERTIAL

Follow the instructions given by your health professional.

This device has been designed for use and application in women of any age and build.

As a general rule, follow the recommendations for use listed below:

- 1) **Empty your bladder**, if necessary.
- 2) **Clean Boltex** with warm water and mild soap.
- 3) **Lubrication:** Use a small amount of lubricating gel in case of vaginal dryness.
- 4) **Insertion:** Relax your pelvic muscles both for inserting and for removing the device.

Bending one or both legs (squatting) can help insertion. Place your index or middle finger on the flat area of the sphere (situated between the removal strings) and push the sphere **into the vaginal passage**. Part of the string will remain outside the body to enable the device to be removed.

5) **Exercising:** Exercise your muscles by walking for an hour or the time specified by your doctor. If you do Kegel exercises as well as using Boltex, the benefits will be greater.

As a general recommendation, start with 15 min/day exercises for 2 weeks. After 2 weeks, you can continue with 30-60 min. exercise sessions.

During treatment, combine fast voluntary contractions (close anus tightly) with slow contractions lasting 15-20 seconds.

Boltex Inertial can be used every day. Improvement in muscle tone is usually noticeable after 6-8 weeks of use. Most women recover their muscle tone after 3-6 months, depending on their degree of initial low muscle tone.

*Note: At first you may experience some muscular discomfort (cramps) from using Boltex Inertial. Should you experience any severe discomfort, see your doctor before*

*continuing use.*



6) **Removal:** Relax your muscles and use the string to pull out the device.

Women suffering from stress urinary incontinence are advised to remove the device in the knee-chest position and with their muscles relaxed, as they have urinary bladder descent and this position makes it easier for them to remove the device in the first few weeks.

Knee-chest position

7) **Cleaning:** Clean the device before and after use with warm water and pH-neutral soap. Store it in its original case or bag to keep it in good condition and prevent contact with dust or other materials.

# Contraindicaciones

## Counter-indications / Contre-indications

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**You should not use Boltex Inertial in the following cases:**

- If you have a vaginal and/or urinary tract infection.
- During menstruation, with absorbent tampons.
- During sex when there is penetration.

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- If you are, or think you may be pregnant, please see your doctor before using the device.
- See your doctor if you feel any pain and/or inflammation after use.
- Do not use for longer than the time recommended by your health professional.
- Carefully follow health & hygiene guidelines before, during and after using the Boltex Inertial.
- Boltex is for personal use only. Do not share.

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